



User Manual
Excel Entice

USER MANUAL

1. INTRODUCTION	2
1.1 IMPORTANT SYMBOLS IN THIS MANUAL.....	2
1.2 USE OF THE SCOOTER	2
2. GENERAL SAFETY NOTES.....	2
2.1 GENERAL SAFETY NOTES.....	2
2.2 SAFETY INFORMATION ON ELECTROMAGNETIC INTERFERENCE	3
2.3 SAFETY INFORMATION ON DRIVING AND FREEWHEEL MODE.....	3
3. THE MOST IMPORTANT PARTS	4
4. DRIVING.....	4
4.1 BEFORE DRIVING FOR THE FIRST TIME	4
4.2 TAKING OBSTACLES	4
4.3 DRIVING UP AND DOWN GRADIENTS.....	5
4.4 PARKING AND STAND STILL	5
5. PUSHING THE SCOOTER BY HAND.....	5
5.1 DISENGAGING MOTORS	6
6. THE CONTROL PANEL	6
6.1 CONTROL PANEL LAYOUT	6
6.2 DRIVING THE SCOOTER	7
6.3 DIAGNOSTICS AND TROUBLE SHOOTING.....	7
7. ADJUSTMENT OPTIONS.....	9
7.1 THE ENTICE CHAIR	9
7.2 ADJUSTING THE SEAT	9
7.3 ADJUSTING THE WIDTH OF THE ARMRESTS	10
7.4 RELEASING THE SEAT SO THAT IT MAY ROTATED AND/OR REMOVED	10
8 ELECTRICAL SYSTEM.....	11
8.1 ELECTRONICS PROTECTION SYSTEM	11
8.2 BATTERIES	11
8.2.2 CHARGING THE BATTERIES.....	12
9. CARE AND MAINTENANCE	13
10. TECHNICAL SPECIFICATIONS.....	14

1. Introduction

Dear user,

Welcome 'aboard' of your new **Entice** scooter!

We wish to thank you for letting us improve your freedom and independence. This model has been designed with your practical needs in mind. It is equipped with modern high tech electronics and special features for a more comfortable ride. Its safety and performance will provide you with years of excellent service and pleasure!

It is from the utmost importance to read this manual very well before using the Excel **Entice** scooter. The manual contains important information concerning the use and small maintenance of the **Entice** . Safe use of the **Entice** is a following issue. This manual should be stored very well. It contains information that might come in handy later!

1.1 Important symbols in this manual



Note: This symbol indicates hints and suggestions which should help make operating the product easier and point out special functions.



Warning: This symbol warns you of danger! Follow the instructions to avoid injury to the user or damage to the product!



Requirements: This symbol indicates a list of the different tools and other requirements you will need to do certain maintenance work.

1.2 Use of the scooter

The intended use of the Excel **Entice** is for people of all ages, who have trouble with walking far distances or walking for a longer time. The **Entice** is an outdoor use scooter. The seat can be adjusted by the user. The seat height can also be adjusted to lower leg length. The seat is able to turn and the flip backward armrests can be used for an easy transfer.

2. General safety notes

2.1 General safety notes



Warning:

Van Os Medical B.V. can not be held responsible for unsafe situations, accidents or damages, if requirements and warnings are not being followed.

- Follow the instructions of this manual
- Do not drive while under the influence of alcohol or medicine.
- Make sure your switch the power system off whenever you get in or out your scooter.

- Be aware that there are only the motor brakes to stop your scooter. When the motors are disengaged, these brakes are automatically deactivated. For this reason, freewheel operation is only recommended on flat surfaces, never on gradients. Never leave your vehicle on a gradient with its motors disengaged. Always re-engage the motors immediately after pushing the vehicle.
- Danger of injury if the On/Off Button is pressed while the vehicle is in motion, due to it coming to an abrupt, sharp stop!
- If you have to brake in an emergency, simply release the drive lever, which will bring you to a halt! Only switch the vehicle off while in motion as a last resort!
- Never sit in the scooter while transported in another vehicle!
- Do not exceed the maximum permissible load (see technical specifications)!
- When maintaining, servicing or lifting any part of your scooter, take into account the weight of the individual components, especially of the batteries! Be sure at all times to adopt the correct lifting posture and ask for assistance if necessary!
- Always use the seatbelt of the scooter.
- Make sure that no injury is incurred by moving parts of the scooter. Make sure that the wheels and other spare parts are fixed.
- Do not connect any electric devices to your vehicle that are not expressly certified by Excel® for this purpose! Have all electrical installations done by your authorized Excel® Dealer!
- Only use original Excel® spare parts, which have been approved for use with this vehicle!

2.2 Safety Information on Electromagnetic Interference

This electric vehicle was successfully tested in accordance with international standards as to its compliance with Electromagnetic Interference (EMI) Regulations. However, electromagnetic fields, such as those generated by radio and television transmitters, and cellular phones, can influence the functions of electric vehicles. Also, the electronics used in our vehicles can generate a low level of electromagnetic interference, which however will remain within the tolerance permitted by law. For these reasons we ask you to please observe the following precautions.

However,

such as those generated by radio and television transmitters, and cellular phones, can influence the functions of electric vehicles.

For these reasons we ask you to please observe the following precautions:

To use the **Entice** close to electromagnetic fields we advise you to switch off the power system. Devices that send electromagnetic fields are telephones, a cellular Phone. Also, the electronics used in our vehicles can generate a low level of electromagnetic interference, which however will remain within the tolerance permitted by law. Is it possible that there will be an unintended event such as a department store alarm going off.

2.3 Safety Information on Driving and Freewheel Mode



Warning:

- Only ever negotiate gradients of up to the maximum defined in the Technical Specifications and only with the backrest in an upright position, and the seat lifter in the lowest position (if installed)!
- Only ever drive downhill at a maximum of 2/3 of the top speed! Avoid abrupt braking or accelerating on gradients!
- If at all possible, avoid driving on slippery surfaces (such as snow, gravel, ice etc.) where there is a danger of you losing control over the vehicle, especially on a gradient! If driving on such a surface is inevitable, then always drive slowly and with the utmost caution!
- Never attempt to overcome an obstacle when on an uphill or downhill gradient!
- Never attempt to drive up or down a flight of steps!
- Always approach obstacles straight on! Ensure that the front wheels and rear wheels move over the obstacle in one stroke, do not stop halfway! Do not exceed the maximum obstacle height (see Technical Specifications)!

- Avoid shifting your centre of gravity as well as abrupt changes of direction when the vehicle is in motion!
- Never use the vehicle to transport more than one person!
- Do not exceed the maximum permissible load!
- When loading the vehicle, always distribute the weight evenly! Always try to keep the centre of gravity of the vehicle in the middle, and as close to the ground as possible!
- Note that the vehicle will brakes or accelerate if you change the Driving Speed while it is in motion!
- Drive through narrow passages in the lowest Driving Speed and with due caution!
-

3. The most important parts

- 1: De-clutching lever
- 2: Release lever for swiveling and removing seat (under the seat on the right)
- 3: Release lever for seat sliding rail adjustment (under the seat on the right in the front).
- 4: Key switch (ON/OFF)
- 5: Primary Charging Socket
- 6: Throttle
- 7: Control panel
- 8: Lever for adjusting the angle of the steering column



4. Driving

4.1 Before driving for the first time

Before you take your first trip, you should familiarize yourself well with the operation of the vehicle and with all operating elements. Take your time to test all functions and driving modes.



Note:

Always use your seatbelt.
When sitting comfortably your drive the safest.

Before each trip, make sure that:

- You are within easy reach of all operating elements.
- The battery charge is sufficient for the distance intended to be covered.
- The seatbelt is in perfect order.

4.2 Taking Obstacles

Your Entice can overcome obstacles and kerbs of up to 8 cm.




Warning:

- Never approach obstacles at an angle!
- Put your backrest into an upright position before climbing an obstacle

Driving up over an obstacle

- Approach the kerb or obstacle slowly head-on. Shortly before the front wheels touch the obstacle, increase the speed and reduce only after the rear wheels have also climbed the obstacle.

Correct

**Driving down off of an obstacle**

- Approach the kerb or obstacle slowly head-on. Before the front wheels touch the obstacle, reduce speed and keep it until also the rear wheels have come down off of the obstacle.

Incorrect

**4.3 Driving up and down gradients**

The Entice can safely climb the following gradients:

- 4-Wheel Version (up to 136 kg payload):12°
- 3-Wheel Version (up to 120 kg payload):12°
- 3-Wheel Version (up to 136 kg payload):10°

 Warning:

- Only ever drive downhill at a maximum of 2/3 of the top speed!
- Always return the backrest of your seat to an upright position before ascending slopes! We recommend that you lean the backrest slightly to the rear before descending slopes!
- Never attempt to ascend or descend a slope on slippery surfaces or where there is a danger of skidding (such as wet pavement, ice etc)!
- Avoid trying to get out of the vehicle on an incline or a gradient!
- Always drive straight in the direction the road or path you are on goes, rather than attempting to zigzag!
- Never attempt to turn around on an incline or a slope!

4.4 Parking and stand still

When parking your vehicle or if your vehicle is stationary for a prolonged period, switch the vehicle's power system off (key switch).

5. Pushing the scooter by hand

The motors of the scooter are equipped with automatic brakes, preventing the scooter from rolling away out of control when the power supply is switched off. When pushing the scooter, the magnetic brakes must be disengaged.

Attention! Switch of the power system Zet het contact af, dit maakt het duwen van de scooter makkelijker.

5.1 Disengaging motors



Warning

Danger of the vehicle running away!

When the motors are disengaged (for push operation), the electromagnetic motor brakes are deactivated! When the vehicle is parked, the levers for engaging and disengaging the motors must without fail be locked firmly into the "DRIVE" position (electromagnetic motor brakes activated)!

The lever for engaging and disengaging the motor is located on the right at the rear of the scooter.

Disengaging the motor

- Make sure the power supply is switched off (key switch).
- Push lever (1) forwards.

Engaging the motor


- Pull lever backwards.



6. The Control panel

6.1 Control Panel layout

1. Seat Lifter activation (if installed)
2. Battery charge indicator
3. Hazard flashers
4. Horn
5. Right turn signal
6. Driving speed adjustment
7. Throttle lever
8. Plug for external charger
9. Reduced Speed Mode
10. Left turn signal
11. Lights



6.1.1 Seat Lifter "Up" and "Down" Buttons

- Press the button to activate the Seat Lifter (if is installed).
- Raise or lower the seat using the throttle lever.



Note:

If the Seat Lifter is not lowered entirely, then the speed of the scooter is automatically reduced to lower the risk of tipping over. This is indicated by the LED just above the buttons being lit. Lower the Lifter to restore speed to it's normal level. The LED will go out.

6.1.2 Status display



Note:

The leftmost diode of the battery charge display serves as an error message display (status display). For an explanation of the Error Codes please see chapter 6.3.

6.1.3 Battery charge display

- All diodes lit: full driving range
- Only red and yellow diodes lit: decreased drive range. Charge batteries at end of journey.
- Only red diodes lit / flashing: battery reserve = very low drive range!
- Charge batteries immediately!



Note:

Total discharge protection: After a certain drive time on reserve battery power the electronics switches the drive off automatically and the scooter will be immobile.

6.2 Driving the Scooter

Switch on the power supply (key switch). The displays on the Control Panel light up. The scooter is ready to drive.



Note:

Note:

If the scooter does not respond after switching on, check the status display (see chapter "Status display").

- Set the driving speed with the speed adjustment knob.
- Gently pull the right driving lever towards you to drive forwards.
- Gently pull the left driving lever towards you to drive backwards.



Note:

The controller is programmed with standard values ex-works. Your Excel® Dealer can program it to fit your requirements.



Note:

To brake quickly, simply let go of the driving lever. It will automatically return to the middle position. The scooter will brake.

6.3 Diagnostics and Trouble Shooting

The electronics system provides diagnostics information to assist technicians to diagnose and correct faults within the scooter system. The existence of a fault will cause the status light to flash in bursts, separated by a pause. The nature of the fault is indicated by the number of flashes in each burst, also referred to as the Flash Code.

Depending on the severity of the fault and impact on user safety, the electronic system will react differently. It may...

- Simply display the Flash Code as a warning and allow normal driving and operation.
- Display the Flash Code, stop the scooter and prevent driving until the electronic system has been turned off and then back on again.
- Display the Flash Code, stop the scooter and prevent driving until the fault has been fixed.

For detailed descriptions of what each Flash Code means, and the probable cause and remedy, see Section 6.3.2 on page 8.

6.3.1 Diagnosing Faults

Use the following troubleshooting guide if the scooter fails to operate.



Note:

Turn the key switch on before beginning any diagnostics.

If the Status Light is OFF
Check that the key switch is turned ON
Check that all cables are connected correctly

If only the leftmost diode of the battery charge display is PERMANENTLY ON
Contact your authorized Excel® Dealer

If the leftmost diode of the battery charge display is FLASHING
Count the number of flashes and refer to the next section.

6.3.2 Flash Codes

Number of flashes	Fault	Impact on Scooter	Notes
1	Battery needs charging	Will Drive	Battery charge is running low. Recharge the batteries as soon as possible.
2	Battery voltage too low	Drive inhibited	Battery charge is empty. Recharge the batteries. If the scooter is left off for a few minutes, battery charge may recover sufficiently to allow driving for a short period of time.
	Lifter raised	Driving speed reduced	Lower lifter completely.
3	Battery voltage too high	Drive inhibited	Battery charge is too high. If a charger is plugged in, unplug it. The electronic system charges the batteries when travelling down slopes or decelerating. Excessive charging in this manner may cause this fault. Turn the scooter power off and then back on again.
4	Current Limit Time Out	Drive inhibited	The scooter has drawn too much current for too long, possibly because the motor has been over-worked, jammed or stalled. Turn the scooter power off, leave it off for a few minutes, and then turn the power back on again. The controller has detected a short-circuited motor. Check the cable harness for short and check the motor: Contact your authorised Excel® Dealer. Check that the declutching lever is in the engaged position.
5	Brake Fault	Drive inhibited	The park brake coil or wiring is faulty. Check the park brake and wiring for open or short circuits. Contact your authorised Excel® Dealer.
6	Out Of Neutral Drive At Power Up	Drive inhibited	Throttle is not in neutral position when turning key switch on. Return throttle to neutral, turn power off, and back on again. Throttle may need to be re-calibrated. Contact your authorised Excel® Dealer.
7	Speed Pot Error	Drive inhibited	The throttle or its wiring may be faulty or incorrectly set up. Contact your authorized Excel® Dealer.
8	Motor Volts Error	Drive inhibited	The motor or its wiring is faulty. Contact your authorised Excel® Dealer.
9	Other Internal Errors	Drive inhibited	Contact your authorized Excel® Dealer.


7. Adjustment options

7.1 The Entice chair

The chair of the Entice is equipped with a number of adjustment possibilities which improves your seat and driving comfort considerably. The entire chair can be height adjusted in a range of 10 cm in 4 steps of 2,5 cm. This adjustment can be done by your dealer.

Adjust backrest angle

The backrest angle is stepless adjustable in a sitting position up to standard 105 degrees . By replacing the hex bolt on the turning point this can be increased up to 135 degrees. The backrest can also be fully flipped forward to ease transport. On the right side of the chair (left shown on picture 36) is a lever. Pull it up and lean forward or rearward to adjust the desired backrest angle . As soon as you reached it you release the lever. The backrest will fixate in your position.

-  **Requirements:**
 Allen key: 4 mm
 Spanner: 10 mm
- Remove the screws on both sides of the seat that hold the backrest (1), using the spanner and Allen key.
 - Adjust the backrest to the desired angle.
 - Reposition the screws and tighten.



7.2 Adjusting the seat

The lever for adjusting the seat back and forth is located under the seat in the front on the right side.

- Pull the lever (1) to release the seat
- Slide the seat forward or backward to the desired position
- Release the lever to lock the seat in place.



7.3 Adjusting the width of the armrests

The knobs to release the armrests are located in the back, under the seat (1).

- Turn the knobs to release the armrests
- Adjust the armrests to the desired width
- Retighten the knobs



7.4 Releasing the seat so that it may rotated and/or removed

The seat can be rotated to the side to ease getting on and off of the scooter. In this position, the seat can also be removed.

The lever for releasing the seat so that it can be rotated is located under the seat on the right (1).

- Push the lever forward to release the seat.
- Rotate the seat to the side.
- If desired, firmly grip the seat by the backrest and the front edge and pull upwards to remove it.



Note:

Do not use the uppermost hole. In this position the seat is too low, and does not have enough clearance above the cowling.

- Adjust the seat height.
- Reposition the bolt and tighten.



8 Electrical system

8.1 Electronics Protection System

The vehicle's electronics are equipped with an overload-protection system.

If the motors are put under considerable strain for a longer period of time (for example, when driving up a steep hill) and especially when the ambient temperature is high, then the electronic system could overheat. In this case the vehicle's power is reduced gradually until it finally comes to a halt. The Status Display shows a corresponding error code. By switching the power supply off and back on again, the error code is cancelled and the electronics are switched back on. It will take approximately five minutes until the electronics have cooled down enough that the motors deliver their full power again.

When the motors are stalled by an insurmountable obstacle, such as a high kerb, and the vehicle driver allows the motors to strain against this hindrance for more than 20 seconds without moving, then the electronics will automatically switch off to prevent the motors from being damaged. The Status Display shows a corresponding error code. By switching off and back on again, the error code is cancelled and the electronics are switched back on.

8.1.1 The main fuse

The entire electric system is protected against overload by two master fuses. The master fuses are mounted on the positive battery cables. They can be replaced only after removing the battery lid.



Note:

A defective main fuse may be replaced only after checking the entire electric system. An Excel® specialised dealer must perform the replacement.

8.2 Batteries

8.2.1 What you need to know about batteries

Power is supplied by two batteries. The batteries are maintenance-free and only need regular charging.

New batteries should always be once fully charged before their first use. New batteries will be at their full capacity after having run through approx. 10 - 20 charging cycles. How fast the batteries will be discharged will depend on many circumstances, such as ambient temperature, condition of the surface of the road, tyre pressure, weight of the driver, way of driving and utilisation of lighting, etc.



Note:

Pay attention to the battery charge indicator! Make sure to charge the batteries when the battery charge indicator shows that battery charge is low. We recommend charging the batteries after each trip, as well as each night over night. Depending on the level of discharge, it can take up to 12 hours until the batteries are fully charged again. Protect your charger for sources of heat such as heaters and direct sunlight. If the battery charger overheats, charging current will be reduced and the charging process delayed.

To avoid damaging the batteries, never allow them to be fully discharged. Do not drive on heavily discharged batteries if it is not absolutely necessary, as this will strain the batteries unduly and greatly shorten their life expectancy.

In case your vehicle is not used for a longer period of time, then the batteries must be charged at least once a month to maintain a full charge. Alternatively, the vehicle can stay connected to the charger. The batteries cannot be overcharged with the specified charger.

Please use only charging devices in Class 2. This class of chargers may be left unattended during charging. All charging devices which are supplied by Excel® comply with these requirements.

8.2.2 Charging the batteries

Make sure you read and understand the battery charger's user's manual, in case one is supplied with it, as well as the safety notes on the front and rear panels of the charger!

Warning:

Danger of explosion and destruction of batteries if the wrong battery charger is used!

- Only ever use the battery charger supplied with your vehicle, or a charger that has been approved by Excel®!

Danger of electric shock and damage to the battery charger if it is allowed to get wet!

- Protect the battery charger from water!

Danger of short circuit and electric shock if the battery charger has been damaged!

- Do not use the battery charger if it has fallen on the ground or been damaged!

Danger of fire and electric shock if a damaged extension cable is used!

- Only ever use an extension cable if it is absolutely necessary! In case you must use one, make sure it is in good condition!

Note:

How can I tell when the batteries are fully charged?

Look through the slot for the de-clutching lever (2) down into the motor and battery compartment.

You will see a green LED. This is the charger LED.

LED ACTIVITY	CHARGING PHASE
LED lights constantly	Start
LED dims and flickers	90% Finished
LED go out / flickers for a few seconds / goes out again	Finished

The Excel® Entice 'Universe' can also be charged using an external charger. The charging socket is on rear edge of the Control Panel (1).

Connecting the charger

- First connect the battery charger to the scooter.
- Then connect the battery charger to the mains power supply and switch on.

Disconnecting the charger

- First disconnect the battery charger from the mains power supply.

Then disconnect the battery charger from the scooter



9. Care and Maintenance



Note:

Have your vehicle checked once a year by an authorised Excel® dealer in order to maintain its driving safety.

Cleaning the vehicle

When cleaning the vehicle, pay attention to the following points:

- Only use a damp cloth and gentle detergent.
- Do not use any scrubbing agents.
- Do not subject the electronic components to any direct contact with water.
- Do not use high-pressure cleaning devices.

Maintenance Jobs	When delivered	Weekly	Monthly
Seat and backrest padding:			
- Check for perfect condition			X
Tyres:			
- Have the tyres checked for specified air pressure (2,5 bar)	X	X	
Front wheels:			
- Front wheels must spin smoothly		X	
- If wheels wobble or do not spin easily, adjust steering pivot pin or front wheel bearing.			X
Rear wheels:			
- Test wheel for firm seat on the axle drive shaft			X
- Rear wheels must spin without wobbling			X
Electronics / Electrical System:			
- Check all plug connections for condition and firm connection			X
- Have batteries been fully charged before the daily operation		Before every trip	
- Are all holders, screws firmly fixed, tight and safe?			X
- Are all electric bulbs of the lightning system (if applicable) in working order?		Before every trip	
Cleaning:			
- Clean all parts carefully		When necessary	

10. Technical specifications

	3 Wheel version	4 Wheel version
Electrical System		
Motor Batteries	295 W / 5100 RPM	295 W / 5100 RPM
Main battery fuse	Standard: 2 x 30 AH Option: 2 x 40 AH 30 A	Standard: 2 x 30 AH Option: 2 x 40 AH 30 A
Charger Input 230V AC, Output 3 A, 24 V DC	On board, Input 230V AC, Output 3 A, 24V DC (external battery charger plug on rear of control panel)	On board, Input 230V AC, Output 3 A, 24V DC (external battery charger plug on rear of control panel)
Weight		
Empty weight (including 30 Ah batteries)	88 kg*	94 kg*
Weight of the heaviest part	21.4 kg	25 kg
Maximum load (payload)	136 kg	136 kg
Dimensions		
Height	120 cm*	120 cm*
Width	63 cm*	63 cm*
Length (w/o bumpers and anti-tippers)	120 cm*	122 cm*
Seat height (measured from the chassis, manually adjustable)	43, 45, 47, 49 cm*	43, 45, 47, 49 cm*
Backrest height (w/o headrest)	43 cm*	43 cm*
Backrest height (/with headrest)	63 tot 69 cm*	63 tot 69 cm*
Backrest angle (manually adjustable)	105°, 120°, 135°, 150°	105°, 120°, 135°, 150°
Seat width	44-48 cm*	44-48 cm*
Seat depth	45 cm*	45 cm*
Armrest height	15-22 cm*	15-22 cm*
Tyres	20-26 cm*	20-26 cm*
Tyre pressure if fitted with pneumatic tyres	3.5 bar (50 psi)	3.5 bar (50 psi)
Driving attributes		
Speed	12 km/h (8 mph)	12 km/h (8 mph)
Safe climbing capability	Tot 136 kg lading: 10° (17%)	Tot 136 kg lading: 12° (21%)
Maximum obstacle height	8 cm	8 cm
Minimum turning radius	120 cm	150 cm
Range according to ISO 7176 **		
With 30 Ah Batteries:	27 km (16.77 miles)	27 km (16.77 miles)
With 40 Ah Batteries:	34 km (21.12 miles)	34 km (21.12 miles)

* Approximately.

** Note: The scooter's range depends strongly on various factors, such battery charge, ambient temperature, local topography, road conditions, tyre pressure, driver's weight, driving habits and the usage the batteries for lighting, actuators, and so on.

CERTIFICATE



EN 12183
EN 1041

ISO 7176-15
ISO 7176

Product identification

Product : *Physically propelled wheelchair*
Brand : *Excel*
Model/type : *Entice*
Version :

Manufacturer:

Name *Van Os Medical B.V.*
Address *Koperslagerij 9*

Country *The Netherlands*

EU Representative:

Name *W. van Os*

Address *Koperslagerij 9*

Country *The Netherlands*

Function *Director*

Technical constructed file

Prepared by: Name *J.M.J. Brouwer BBA*
Function: *Research and Development*
Issue date: *01 - 12 - 2010*
TCF date: *01 - 12 - 2010*
Recertification date: -

Certificate/report no.: N/A

Crash test report no. according to ISO 7176/19: N/A

Means of conformity

The product is in conformity with Directive 93/42/EEG based on the use of a Technical construction file in accordance with Article 9 (Class I products) of the Directive

Signature of EU representative:

Place : *Steenbergen*

Date : *01 - 12 - 2010*

Number: *VOS.TCF.EX.0906*



Van Os Medical UK Ltd.

Excel House
Escrick Business Park
Escrick, York
YO19 6FD

Tel. +44 (0) 1904 720170
Fax. +44 (0) 1904 720398

E-mail: sales@vanosmedical.com
www.vanosmedical.com