

Jazzy

POWER CHAIRS

Jazzy

600 ES

A unique power chair, with exceptional comfort and support.



Images illustrated:

1. Fold Flat Seat
2. Active-Trac® ATX Suspension
3. Exclusive Mid-Wheel Drive design
4. Width adjustable arm rests and Comfort High Back Seat
5. J600ES with Comfort High Back Seat

* Images used for illustration purposes only



CRASH TESTED



Introduction

Offering unique comfort, support and heightened suspension. The Jazzy 600 ES is an exceptional indoor/outdoor powered wheel chair. Offering superb, all-around performance and is a good power chair choice.

The Jazzy 600 ES features a great selection of standard features. Such as a comfort high back seat, ATX Suspension, along with smooth all-round performance from the mid-wheel drive system. Enabling excellent durability and outstanding functionality.

Standard Features



Mid-Wheel Drive design, which is ideal for tight turning circles



Luxurious comfort high back seat



ATX suspension



Sophisticated modern design

Technical Specifications

Drive Wheels (cm) 35.56 cm solid tyres only	OMNI Casters (cm) Front and Rear: 15.25 cm	Max Speed (cm) 4 Mph (6.4 km/h)	Ground Clearance (cm) 7 cm At battery tray	Turning Radius (cm) 52.1 cm	Overall Length (cm) 102.2 cm	Overall Width (cm) 62.2 cm	Seating Options (inches) 16"W-AcJ 16" - 18"D 18"W-AcJ 18" - 20"D 20"W-AcJ 18" - 20"D	Batteries x 2 volt 50 - 55 AH Deep Cycle	Available Electronics 70A PG VR2 4 or 6 key	Weight Capacity (stone) 21 Stone 6 lbs (137kg)
Base Weight (kg) 56.1 kg	Battery Weight (kg) 16.9 kg	Colours: Black Blue Red		Warranty Lifetime limited warranty on frame; 1-year limited warranty on electronics; 18-month warranty on drive control system; 2-year limited warranty on drivetrain						



Tel: 01869 324600
 Email: sales@quantumrehab.co.uk
 Web: www.quantumrehab.co.uk
 (Version 2: August 2015)

Notes

- Speed varies with user weight, terrain type, battery charge, battery condition and tyre pressure.
- Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.
- Battery weight may vary based on manufacturer.