



Salsa® R



Salsa® R

Compact yet powerful...
... at home outside as well as inside!



Great outdoor performance

- Up to 10cm (4") kerb climbing
- Speed 4mph (6kph) (6mph (10kph) optional)
- 60 Ah Gel batteries provide up to 20 miles (32km) range
- 8° dynamic stability (14%)
- Lights (optional)
- 14" drive wheels
- 9" castor wheels

Excellent indoor manoeuvrability

- Compact turning radius of 85cm (33.5")
- 58cm (22.8") width with 12.5" wheel option gives easy access through doors and confined spaces
- Low seat height for easy table and vehicle access (lowest 37cm (14.6") with 0° manual tilt)

Fully adjustable seating system with retro-fit options for easy assessment, refurbishment and servicing

One fully adjustable seat with more than TWO MILLION seating combinations as standard!

- Only four tools required for adjusting the complete seating system
- Five seat widths and five seat depths
- Five seat heights and four seat angles
- Three hanger depths plus adjustable footplates
- Full height, depth and angle adjustment on flip-back armrests
- Six backrest angles and three backrest heights

Safety on slopes, ramps and in transport

- Fully transit tested (ISO 7176 Part 19) with six point tie downs
- Quick release, removable backrest for easy transport
- Centre of Gravity seat adjustment for optimised stability
- Tested for use with 23cm (9") dual ramps



NEW: Low push force, emergency freewheel allows the Salsa to be moved very easily when it is in freewheel



compatible



Flip-up armrests and footrests for easy accessibility



Lift-off backrest for easy transportation and storage

Optional Pro Seating: Comfort Seating System

JAY COMFORT

A variety of Comfort Seating options are available, these include:

- Comfort seat pan
- Standard and optimum comfort cushion
- Shallow and deep contour backrests
- Comfort calf pads
- Comfort and flexible headrests



Optimum comfort cushion



Deep contour backrest



Comfort headrests

Advanced electronics

VR-2 electronics as standard

- 70 Ah controller as standard. 90 Ah heavy duty version optional
- Dual control (option)

R-net expandable controls (optional)

- Provides simple programming options
- Patented assignable buttons functionality
- Controls up to six actuator controls through the joystick
- Allows unique programming of profiles
- Bluetooth mouse control
- Dual control

Full range of powered options available

Available powered options:

- C.O.G tilt-in-space
- Lift and tilt
- Recline
- ELR/ALR powered legrests
- Easy adaptation of powered options if required (upgrade kit)

Optional retro-fit kits:

A full range of upgrade kits can be added to your chair at a later date if your needs change

Salsa®

Freedom indoor and out



Technical Data: Seating

Seat depth:	41, 43, 46, 48, 51cm (adj 41 - 51cm)	16, 17, 18, 19, 20"
Seat width:	41, 43, 46, 48, 51cm (adj 41 - 51cm)	16, 17, 18, 19, 20"
Seat height:	37 - 46cm (no tilt), 41 - 44cm (with tilt), 44 - 44.5cm (lift & tilt)	14.5 - 18" (no tilt), 16 - 17.5" (with tilt), 17.5 - 18" (lift & tilt)
Backrest angle:	-4°, 0°, 4°, 8°, 12°, 16°	
Backrest width:	41, 43, 46, 48, 51cm (adj 41 - 51cm)	16, 17, 18, 19, 20"
Lower leg length:	33 - 51cm	13 - 20"
Legrest angle:	70°, ALR, ELR	
Armrest height:	Adj. from 22 to 31cm	8.5 to 12"
Armrest depth:	Adj. from 31 to 41cm	12 to 16"

Specification

Overall width*:	58, 60, 62cm	23, 24, 24.5"
Overall length:	114cm	45"
Chair weight including batteries**:	113 - 165kg	250 - 363lbs
Turning radius:	85cm (12" wheels)	33" (12" wheels)
Battery capacity	60Ah gel batteries	
Range:	32km	20 miles
Speed:	6kph, 10kph	4mph, 6mph
Kerb climbing ability***:	8.5cm (8" castors), 10cm (9" castors)	3.3" (8" castors), 4" (9" castors)
Maximum user weight:	140kg	22 st
Max safe slope:	8°	14%
Transit tested to ISO 7176 Part 19 and passed:	Yes	
Warranty:	2 years on frame and electronics, 3 years on motors, 1 year on batteries (conditions apply) and chargers	

Find out more today!

Visit www.quickiepower.com for more information, to find your local dealer or to request a demonstration.

For further information on the full specification, options and accessories, please refer to the order form. All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have. If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.com or alternatively is available on request in large text.



Get SOCIAL with QUICKIE

Keep up-to-date with the latest news on the Quickie Salsa R and other Quickie products through our Facebook page:

www.facebook.com/QuickieWheelchairsUK



Sunrise Medical Limited
Thorns Road
Brierley Hill
West Midlands
DY5 2BQ
Tel: +44 (0) 845 605 66 88
Fax: +44 (0) 845 605 66 89
www.SunriseMedical.com

* 12" drive wheels, 14" drive wheels, low push to force drive wheels

** Please note: Total weight of the powerchair depends on the configuration

*** With kerb climber