

# Hybrid-Turn

The Hybrid-Turn mattress is a unique hybrid concept. This mattress offers the traditional benefits of high specification visco-elastic foam and true alternating therapy by simply rotating the mattress. Suitable for clients weighing up to 254kg / 40st who are at high risk of developing pressure ulcers.

In one easy turn, this mattress allows for immediate step down or step up to a full body length alternating surface.



Product Specification			
Product Code	1001280001	Pump Dimensions	254 x 130 x 95mm
Max. Weight	254kg / 40st	Pump Weight	1.6kg
Risk Factor	Up to High Risk	Crib Rating	Source 0, 1 & Crib 5
Operating Cycle	10, 15, 20, 25 mins	Power Requirements	
Mattress Dimensions (Inflated)	1980 x 880 x 160mm	Safety Standards	BS5852 & BS7177
Foam & Cover Warranty	5 years	Voltage	AC230V 50 Hz 0.5A
Pump Warranty	2 years	Noise Level	≤30db
Mattress Weight	20kg	Fuse Rating	T1AL/250v



# The Mattress

The Hybrid-Turn mattress combines both static and alternating therapy, thus only one support surface is required.



Alternating side

## Alternating Side

- Alternating cells relieve the pressure and assist with pressure ulcer management
- Suitable for up to Category 4 pressure ulcers

**Step up** – Alternating surface (pressure relieving)



Static side

## Static Side

- High specification visco-elastic foam conforms to the body, reducing the pressure for effective pressure ulcer prevention
- Suitable for up to Category 2 pressure ulcers

**Step down** – Static surface (pressure reducing)

## Sodium Hypochlorite Resistant Cover

- Can withstand a higher concentration of cleaning substances
- Offers greater resistance to abrasion
- High frequency welded seams protect against fluid ingress
- 360° zip to prolong longevity
- White inner substrate to aid infection control guidelines

## Hose Management

Integral covered hoses are stored in convenient storage hoops



## CPR

The easy access and simple to use CPR facility tag allows for the rapid deflation of the mattress.



## Individually removable cells

The cell construction enables modular components to be individually replaced or cleaned, thus helping to reduce ongoing costs.



### Registered Office:

Park House Healthcare Ltd, Whitehall 26 Business Park, 1 Heathfield Lane, Birkenshaw, West Yorkshire, BD11 2HW  
Company Reg. No. 1831906 V.A.T. Reg. No 406 0032 16

